


# Panasonic®

## Operating Instructions Software Guide

### TLD READER

Model No. UD-7900N14EU



Before attempting to connect or operate this product, please read these instructions carefully and save this manual for future use.

The model number is abbreviated in some descriptions in this manual.

Refer to the manuals provided with the devices in use or descriptions of how to perform the basic operations relating to the Windows® operating systems and how to handle the PC and peripheral devices.

## Trademarks and Registered Trademarks

- Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.
- Intel and Core i5 are trademarks or registered trademarks of Intel Corporation or its subsidiaries in the United States and other countries.
- IBM and PC/AT are trademarks of International Business Machines Corporation, registered in many jurisdictions worldwide.
- Macintosh is a registered trademark of Apple Inc.
- All other trademarks identified herein are the property of their respective owners.

## Prohibition/Limitation of Liability

- IN NO EVENT SHALL Panasonic System Networks Co., Ltd. BE LIABLE TO ANY PARTY OR ANY PERSON, FOR THE CASES, INCLUDING BUT NOT LIMITED TO ANY ACCIDENT OR DAMAGE CAUSED BY INAPPROPRIATE USE, NEGLIGENT OPERATION OR MISTAKEN SETUP OF THE USER
- IN NO EVENT SHALL Panasonic System Networks Co., Ltd. BE LIABLE TO ANY PARTY OR ANY PERSON, FOR THE CASES, INCLUDING BUT NOT LIMITED TO ANY ACCIDENT OR DAMAGE CAUSED BY DISASSEMBLE, REPAIR OR MODIFICATION OF THE PRODUCT BY A PERSON OR ENTITY OTHER THAN Panasonic System Networks Co., Ltd. OR THE SERVICE COMPANY DESIGNATED BY Panasonic System Networks Co., Ltd.
- No part of this book may be reproduced or used without the prior permission.
- Some of the descriptions in this document may differ from the data at the time of purchase since they are based on the data at the time of documentation.
- Contact Panasonic System Networks Co., Ltd. for questions.

## Precautions

Be sure to refer to the precautions described in the manuals of the PC to install this software.

## Safety Instructions

Be sure to refer to the safety instructions described in the manuals of the PC to install this software.

# Contents

Trademarks and Registered Trademarks.....	1
Prohibition/Limitation of Liability .....	1
Precautions .....	1
Safety Instructions.....	1
Contents .....	2
1. Operation .....	4
■ Launching of application software.....	4
■ How to start measurement.....	5
■ Create a batch file .....	5
■ Set the magazine .....	6
■ Status and display .....	7
2. Reviewing Data.....	9
Exit button .....	9
■ Batch file to be applied.....	9
■ Review data .....	10
■ Export Data.....	11
■ Graph output.....	15
Description of Menus .....	17
■ Main screen .....	17
■ Menu Bar .....	20
■ Menu Bar File Command.....	21
New/Batch.....	21
Exit UD-7900.....	21
■ Menu Bar View Commands.....	22
Badge Results .....	22
Status Results .....	23
Reports.....	23
Re-Output Measured Data.....	23
■ Menu Bar Mode Commands.....	24
Reader Check/Loader .....	24
Reader Check/Changer .....	24
Reader Check/Lamp F3 and F4 Check .....	25
■ Menubar Parameters Commands.....	26
List Parameters .....	26
Processing Modes .....	27
Timing Parameters .....	29
Reader Functions (710 Read Mode) .....	30
Reader Functions (716 Read Mode) .....	31
Calibration Settings.....	32
Heat Lamp Power.....	33
RS232 Options .....	33
Parameter Dump to Computer.....	34
Parameter Dump to File.....	34
■ Menu Bar Help.....	35
About UD-7900.....	35
3. How to Operate Windows .....	36
■ Menu bar navigation.....	36
Appendix-A Error Code.....	37
■ Error code table .....	37

Appendix-B Example of Report.....39  
    Export Data.....39  
    Badge Results .....39  
    Tray Summary.....40  
    Status.....40  
Appendix-C Communication Specifications .....41  
    Basic specifications .....41  
    710/791 Interface mode .....42  
    716 Interface mode.....44  
Appendix-D Parameter list.....46

# 1. Operation

The TLD reader UD-7900P automatically measures a total of 500 badges (10 magazines, 50 badges per magazine) inserted in magazines.

Measurement data will be stored in a file and can easily be checked, analyzed and used to create reports.

## ■ Launching of application software

After checking the setup of the system component devices (connection of signal cables and turning on the power), launch the application software.

When the application software is launched, the main screen (Fig. 1-1) will be displayed. Processing will be executed continuously from the previous badge. The current badge position will be displayed in the current badge information (Fig. 1-2).

If there is no batch at the application software launch, “DefaultBatch” will automatically be created.




Fig. 1-1 Main screen





Fig. 1-2 Current badge information

In the “Location” field, the magazine number and the badge number will be displayed.  
(Example): Values 01-00 in the “Location” field indicate that magazine number “1” (shelf position 1) is waiting for being loaded.

## ■ How to start measurement



Measurement of all badges will start by clicking the [START] button on the main screen. Measurement will stop when the [STOP] button is clicked.

Measurement will also automatically stop when an error occurred. In this case, removing of the error cause and restarting of the measurement will make possible to write the corresponding data.

By clicking the [Reset Motors] button, the mechanical initialization operation will start.

Magazines located in the shelf of the magazine changer will automatically be loaded in order during the badge measurement.

This product will skip shelves without magazines and shift to the next shelf. When processing of the last shelf is completed, measurement will be ended and the changer shelf will be moved to the initial position.

When the magazine with badges is taken in, the badge of the top slot will be sent to the number read position and ID code, badge code and rank code will be read. Then, the dark count and other data will be read and measurement of each element will be executed.

When the measurement is completed, the badge will be returned to the magazine from which the badge originally came and the magazine will be advanced by 1 slot.

When measurement of 1 slot is completed, the magazine will be returned to the shelf and the shelf will be advanced by 1 step. After this, these processes will be repeated.

If a mechanical error has occurred, the measurement cannot be started. Therefore, press the [Reset Motors] button to initialize.

## ■ Create a batch file

First, create a batch file to manage all measurement data at once (It is unnecessary to create it if continuously use the old batch file or only check the file or data.). It is crucial to newly create a batch file to prevent the old batch file data from being lost.

To create a batch file, select [File] - [New] - [Batch] from the menu bar. On the “New Badge” screen that is displayed at this time, register the badge ID, operator and comment (Fig. 1-3 and Fig. 1-4). By clicking the [Make Batch] button, a new file will be created and the screen will return to the main screen. The information registered at this time will be displayed on the “Status Log” window (Fig. 1-5).




Fig. 1-3 New batch




Fig. 1-4 New batch (example)




Fig. 1-5 Program status


## ■ Set the magazine

When executing a measurement, firstly set the magazine to the automatic changer. Refer to the operating instructions of UD-7900 for how to set the magazine.

If there is an unnecessary magazine in the measurement section, remove it by clicking [Unload Tray]. If the shelf position is not “1”, return the automatic changer to shelf position 1 by clicking [Changer Rev.].

The current badge information will display the shelf position in the form of “[Shelf number] - [Slot number]”. If the shelf number is “1” and there is no magazine in the reader, “Batch Status” will display 01-00.

If the magazine is not in the measurement section, “00” will be displayed for the slot position.



Changer Rev	Reversely advances the shelf by 1 step.
Changer Adv	Advances the shelf by 1 step.
Mag Adv	Advances the magazine by 1 slot.
Mag Rev	Reversely advances the magazine by 1 slot.
Unload Tray	Returns the magazine from the reader to the automatic changer.
Load Tray	Takes the magazine out of the automatic changer and send it to the reader.
START	Starts the measurement of a badge.
STOP	Stops the measurement of a badge.
Reset Motors	Executes the mechanical initialization operation.

Magazines at the bottom shelf will be set first. The bottom stage of the shelf is shelf position 1 and the top stage of the shelf is shelf position 10.

## ■ Status and display

Statuses and data during measurement will be displayed. The above includes each status of device, program, power source, badge and operation, analog input data (option) and glow data.

## Measured Data

Each reading value and calculation result during the badge measurement will be displayed. The “Batch” display and the “Location” display are useful to check the magazine position when manually controlling the automatic changer.

## Status Log

Log information including start of measurement, end of measurement, unmatched status, etc. will be displayed. It is possible to browse all log information even after completing the measurement. In this case, use the scroll bar on the right side of the window. The log information will be stored in the “batch” directory for reference in future.

The latest error message will be displayed in red.



## Lamp Status

---

Voltage					
Lamp V.	<input type="text"/>	V	Lamp B.V.	<input type="text"/>	V
Lamp C.	<input type="text"/>	A	Lamp B.C.	<input type="text"/>	A
Lamp P.	<input type="text"/>	W	Lamp B.P.	<input type="text"/>	W


Lamp Flash Count
<input type="text"/>

Each voltage value will be displayed.

The lamp voltage, current and power will be displayed when the measurement is being executed or when "Reader Check" or "Lamp F3/F4 Check" is being executed.

## Glow Curve Display


---



Glow data at measurement will be displayed in auto scale.

## 2. Reviewing Data

The following are descriptions of how to check data and of how to create each report. It is also possible to create a report from the latest batch file or a batch file used in the past. Select [View/Reports] on the main screen.



**Exit button**

**View/Reports – Exit Command**

Close the data review screen and return to the main screen.

### ■ Batch file to be applied

At the default, the latest batch file is selected to be applied (Fig. 2-1). To change to a batch file used in the past, select a batch ID from the list (Fig. 2-2). After selecting the batch file, update the display to the content of the corresponding batch file.




Fig. 2-1 Select the latest batch ID file




Fig. 2-2 Select from the Batch ID list

## ■ Review data

After the measurement, it is possible to review data in text format or in graph format .

### Text format

**Export Data**  
(Refer to pages 12 and 40.)

“Export Data” creates detailed reports of each element.

**Badge Results**  
(Refer to pages 13 and 40.)

“Badge Results” displays all measurement results in the batch.

**Tray Summary**  
(Refer to pages 14 and 41.)


“Tray Summary” displays the average value and calculation result of “o” other than measurement values.

**Status Log**  
(Refer to pages 15 )

“Status Log” displays errors occurred in the selected batch.


### Graph format

**Badge Information**  
(Refer to page 16.)



“Badge Information” displays glow curves of the selected badge.

**Element Information**  
(Refer to page 17.)



“Element Information” displays related parameters other than the P counter information of the specified element.

## ■ Export Data

“Export Data” creates a detailed report of 4 elements on all measurements. This report can be imported for spread sheet software and also be cut & paste on word processing software. To create this report, select [File≠Export Data] from the menu of the “Review Data” screen.

### Format of “Export data”

DateTime				Badge measurement time
Try				Magazine number
Slit				Slot number
BdgNo				Badge number
Typ				Badge type
Err				Badge error code
E1Dose	E2Dose	E3Dose	E4Dose	Element # Dose value
E1Cnts	E2Cnts	E3Cnts	E4Cnts	Element # Count
E1Pre	E2Pre	E3Pre	E4Pre	Element # Pre
E1Main	E2Main	E3Main	E4Main	Element # Main
E1Post	E2Post	E3Post	E4Post	Element # Post
E1PreCnts	E2PreCnts	E3PreCnts	E4PreCnts	Element # Pre-count
E1MainCnts	E2MainCnts	E3MainCnts	E4MainCnts	Element # Main count
E1PostCnts	E2PostCnts	E3PostCnts	E4PostCnts	Element # Post-count
E1Err	E2Err	E3Err	E4Err	Element # Error code
				Drk
				Dark count
				Ref.Element
				Reference element
				Ref.Light
				Reference light
				C
				Rank C
				E
				Rank E
				Unit
				Unit of dose

## Badge Results

“Badge Results” provides the following information about the selected badge. To create this report, select [View≠Badge Results] from the menu of the “Review Data” screen.

Item	Description
Batch:	Batch ID
Operational Time:	Date and time when measurement started
Operator:	Batch operator
Comment:	Comment
Time	Badge measurement time
Try	Magazine number
Slit	Slot number
BdgNo	Badge number
Typ	Badge type
Err	Badge error code
E1	Element 1 dose value
E2	Element 2 dose value
E3	Element 3 dose value
E4	Element 4 dose value
Drk	Dark count
Ref	Reference element
Corr	Reader sensitivity correction factor
C	Rank C
E	Rank E

### Example of output

```


Batch DefaultBatch Processing Report
Operational Time: 12/11/2016 3:46:06 PM Operator:
Comment:

Time Try Slit BadgeNo Typ Err E1 E2 E3 E4 Drk Ref Corr C E
15:46:06 01 03 0702345 10 18uSv 18uSv 185uSv 8uSv 002 017 0.450 7 2
15:46:36 01 04 6101379 10 18uSv 22uSv 221uSv 9uSv 003 000 0.450 5 1
15:46:56 01 05 6101338 10 29uSv 17uSv 155uSv 7uSv 004 000 0.450 2 1
15:47:16 01 06 6101042 10 41uSv 30uSv 209uSv 11uSv 001 012 0.450 0 0
15:47:35 01 07 0702138 10 22uSv 14uSv 140uSv 7uSv 005 000 0.450 5 1
15:47:54 01 08 0702466 10 28uSv 28uSv 278uSv 9uSv 003 010 0.450 0 3
15:48:14 01 09 0702272 10 30uSv 62uSv 260uSv 9uSv 002 007 0.450 0 2
15:48:33 01 10 0702410 10 60uSv 35uSv 208uSv 2uSv 001 009 0.450 7 7

```

# Tray Summary

In “Tray Summary”, it is possible to check the element dose value, average value and “o” (by magazine) of the selected batch. To display this report, select [Tray Summary/Tray #] (#: 1-10) on the “Review Data” screen.



## Example of output

Batch DefaultBatch Tray Report  
Operational Time: 12/11/2016 4:02:59 PM Operator:  
Comment:

Sl#	BadgeNo	E1	E2	E3	E4	Average	Sigma
01	6100003	34uSv	16uSv	125uSv	8uSv	45uSv	46.71
02	6101248	13uSv	11uSv	165uSv	7uSv	49uSv	67.01
03	0702345	18uSv	18uSv	185uSv	8uSv	57uSv	73.87
04	6101379	18uSv	22uSv	221uSv	9uSv	67uSv	88.75
05	6101338	29uSv	17uSv	155uSv	7uSv	52uSv	59.97
06	6101042	41uSv	30uSv	209uSv	11uSv	72uSv	79.39
07	0702138	22uSv	14uSv	140uSv	7uSv	45uSv	54.67
08	0702466	28uSv	28uSv	278uSv	9uSv	85uSv	111.27
09	0702272	30uSv	62uSv	260uSv	9uSv	90uSv	99.81
10	0702410	60uSv	35uSv	208uSv	2uSv	76uSv	78.80
...	...	...	...	...	...	...	...
41	0702157	34uSv	47uSv	177uSv	17uSv	68uSv	63.40
42	0702340	60uSv	22uSv	252uSv	15uSv	87uSv	96.65
43	6100198	24uSv	13uSv	305uSv	18uSv	90uSv	124.19
44	0702415	80uSv	28uSv	301uSv	20uSv	107uSv	114.21
45	6101267	36uSv	55uSv	273uSv	20uSv	96uSv	102.94
46	6100999	34uSv	34uSv	180uSv	15uSv	65uSv	66.42
47	6101063	61uSv	15uSv	263uSv	16uSv	88uSv	102.30
48	6301820	32uSv	8uSv	178uSv	15uSv	58uSv	69.69
49	0702215	29uSv	35uSv	292uSv	14uSv	92uSv	115.44
50	0203013	0uSv	0uSv	0uSv	0uSv	0uSv	0.00
No. of Data		050	050	050	050		
Average		32uSv	30uSv	188uSv	10uSv		
Sigma		16.65	24.01	67.08	4.60		

Item	Description
Batch:	Batch ID
Tray:	Magazine number
Operational Time:	Date and time when measurement started
Operator:	Batch operator
Slit	Slot number
Badge No	Badge number
Element (1-4)	Element measurement value
(in line) Average	Measurement average value (by batch)
(in line) Sigma	$\sigma$ (by batch)
No. of Data	Number of element
(bottom) Average	Measurement average value (by element)
(bottom) Sigma	$\sigma$ (by element)

## Status Log

---

In “Status Log”, it is possible to check errors occurred in the selected batch. To display this report, select [View/Status Log] on the “Review data” screen.

### Example of output

#### TLDErrorLog.txt

```

2016/12/11 15:30:47 UD-7900 Program Started with 710 Read Mode.
2016/12/11 15:30:49 Batch DefaultBatch created.
2016/12/11 15:32:18 START button pressed.
2016/12/11 15:32:20 Processing started on batch DefaultBatch: tray 1.
2016/12/11 15:32:58 STOP button pressed.
2016/12/11 15:37:33 UD-7900 Program Shutdown.


```

## ■ Graph output

In “Graph Output”, mainly glow curves and the related parameters will be displayed.

## Badge Information

Badge Information: Displays glow curves and the related parameters.



In “Badge Information”, glow curve data can be displayed. Select TRAY of badge to be displayed and BADGE.



## Element Information

Element Information: Displays glow curves and the related parameters.

View Reports Window

File Tray Summary

Batch Information

Batch ID: DefaultBatch

Operator: \_\_\_\_\_

Date: 2016/12/11 15:30:47

Comment: \_\_\_\_\_

Badge

Location: 1 33

ID: 0702889

Type: UD810

Rank E-C: 1 - 1

Error: \_\_\_\_\_

Count

Dark: 0001

Ref.Element: 0019

Ref.Light: 2122

Exit

---

Badge Information | Element Information

Element No.: 1

Element Type: 7Li11BO

Dose

Dose: 14uSv

Pre: 7uSv

Main: 14uSv

Post: 467uSv

Count

Pre: 0133


Main: 0265

Post: 1902

Save as

Wave Form

Photon Counter



In “Element Information”, glow curve data can be displayed. Select Element No. of badge to be displayed and BADGE.

# Description of Menus

Descriptions about details of screens and operations are provided here. First, descriptions about content to be displayed on the main screen, and then about each menu are provided. These contain the status displays and the controls to operate through these displays.

## Main screen

Badge Read Window

File View Mode Parameters Help

Status Log

3/9/2017 4:18:47 PM Processing started on batch 030904: tray 4.  
 3/9/2017 4:34:26 PM Processing started on batch 030904: tray 5.  
 3/9/2017 4:50:12 PM Processing started on batch 030904: tray 6.  
 3/9/2017 5:05:53 PM Processing started on batch 030904: tray 7.  
 3/9/2017 5:21:30 PM Processing started on batch 030904: tray 8.  
 3/9/2017 5:37:57 PM Processing started on batch 030904: tray 9.  
 3/9/2017 5:53:42 PM Processing started on batch 030904: tray 10.  
 3/9/2017 6:09:44 PM Processing stopped on batch 030904.

Batch ID: 030904 Location: 1 - 0 Cor. Factor: 0.424 Lamp Flash Count: 8296 Door:  Open

Location	Badge ID	Type	E/C	Dark Cnt	Ref. Ele.	Ref. Light	Element1	Element2	Element3	Element4
1-1	0000194	800	0/0	8	231	0.421	61uSv	64uSv	2uSv	2uSv
1-2	0000195	800	0/0	9	173	0.424	71uSv	14uSv	2uSv	3uSv
1-3	0000198	800	0/0	14	120	0.418	7uSv	32uSv	1uSv	1uSv
1-4	0000521	800	0/0	12	145	0.423	35uSv	32uSv	2uSv	1uSv
1-5	0000152	800	0/0	13	125	0.425	42uSv	4uSv	2uSv	1uSv
1-6	0000200	800	0/0	11	155	0.419	47uSv	18uSv	0uSv	2uSv
1-7	0000199	800	0/0	6	174	0.423	32uSv	39uSv	2uSv	2uSv
1-8	0000190	800	0/0	14	130	0.424	21uSv	1.63mSv	0uSv	0uSv
1-9	0000191	800	0/0	14	122	0.416	4uSv	18uSv	1uSv	0uSv
1-10	0000192	800	0/0	15	146	0.423	21uSv	21uSv	1uSv	2uSv
1-11	0000193	800	0/0	9	118	0.421	53uSv	36uSv	1uSv	2uSv

Glow Curve

Element1 Element2 Element3 Element4

Ready

Voltage  
 Lamp V.: 14.38 V Lamp B.V.: 01.26 V  
 Lamp C.: 20.20 A Lamp B.C.: 02.44 A  
 Lamp P.: 290 W Lamp B.P.: 3 W

Mode  
 NCR Parity  Anneal Mode

Operation  
 Changer Rev.  
 Mag. Adv. Mag. Rev.  
 Changer Adv.  
 Load Tray Unload Tray  
 START STOP  
 Reset Motors

The main screen is the initial screen when operating UD-7900. With the main screen as the starting-point, continuous work with the current batch, creation of new batch, data review of old batch, data analysis, creation of report, etc. can be executed.

Status Log

3/9/2017 3:17:47 PM START button pressed.  
 3/9/2017 3:17:49 PM Processing started on batch 030903: tray 1.  
 3/9/2017 3:30:32 PM STOP button pressed.  
 3/9/2017 3:30:42 PM Processing stopped on batch 030903.  
 3/9/2017 3:31:16 PM Batch 030904 created.  
 3/9/2017 3:31:32 PM START button pressed.  
 3/9/2017 3:31:34 PM Processing started on batch 030904: tray 1.  
 3/9/2017 3:34:05 PM D5 Control limit post dose (P.93) exceeded on element 2 at location 1-8.

“Status Log” displays start of measurement, end of measurement and error/abnormal data occurrence and their occurrence time and date. It is possible to display the past content (However, this is available only after launching the application software. When the application software is stopped once and restarted, content will be initialized.) using the scroll bar on the right side of the screen.

## Lamp Status

Voltage Lamp V. <input type="text"/> V    Lamp B.V. <input type="text"/> V Lamp C. <input type="text"/> A    Lamp B.C. <input type="text"/> A Lamp P. <input type="text"/> W    Lamp B.P. <input type="text"/> W	Lamp Voltage¥B.V.	Lamp voltage/Bleeder voltage (average)
	Lamp Current¥B.C.	Lamp current/Bleeder current (average)
	Lamp Power¥B.P.	Lamp power/Bleeder power (average)
	Lamp Flash Count	Cumulative lamp flash count
Lamp Flash Count <input type="text"/>		

“Lamp Status” displays the lamp voltage, lamp current, lamp power and lamp flash count.

## Batch ID

Batch ID	Location	Cor. Factor
<input type="text" value="new batch"/>	<input type="text" value="1"/> - <input type="text" value="11"/>	<input type="text" value="12.4"/>

Header	Description
Batch ID	Batch ID
Location	Current magazine number and slot number
Cor.Factor	(10 times CAL) count

The current badge ID information will be displayed for “Badge ID”. Displaying of measurement result by each item will go on even when the badge measurement processing is in progress. “Location” displays the magazine position even when the magazine is being sent manually.


## Measured Data

Location	Badge ID	Type	E/C	Dark Cnt	Ref. Ele.	Ref.Light	Element 1	Element2	Element3	Element4
1-1	0000000	802	7/1	3	26	0.441	0uSv	0uSv	0uSv	0uSv
1-2	0320009	802	0/0	1	22	0.441	5uSv	5uSv	0uSv	0uSv
1-3	0000000	802	6/0	2	5	0.441	3uSv	14uSv	0uSv	0uSv
1-1	0000000	800	0/0	1	16	0.444	8uSv	8uSv	1uSv	0uSv
1-2	0000000	800	0/0	2	6	0.444	0uSv	1uSv	0uSv	0uSv
1-3	0000000	800	0/0	4	9	0.441	33uSv	29uSv	43uSv	44uSv
1-4	0000000	800	0/0	0	19	0.441	6uSv	12uSv	1uSv	1uSv
1-5	0000000	800	0/0	2	5	0.440	8uSv	8uSv	0uSv	0uSv
1-6	0000000	800	0/0	3	4	0.440	0uSv	0uSv	0uSv	1uSv
1-7	0000000	800	0/0	0	31	0.435	6uSv	6uSv	0uSv	1uSv
1-8	0000000	800	0/0	1	21	0.435	20uSv	20uSv	43uSv	41uSv

Header	Description
Location	Current magazine number and slot number
Badge ID	TLD Badge number
Type	TLD Badge type
E/C	Rank E value/Rank C value
Dark Cnt	Dark count
Ref.Ele	Reference element
Ref.Light	Reference light
Element 1	Element 1 measurement result
Element 2	Element 2 measurement result
Element 3	Element 3 measurement result
Element 4	Element 4 measurement result


In “Measured Data”, displaying of measurement result by each item will go on even when the badge measurement processing is in progress.

## Glow Curve Display




“Glow Curve Display” displays preview of glow data.

## Autochanger Control



Changer Rev	Reversely advances the shelf of the automatic changer by 1 step.
Changer Adv	Advances the shelf of the automatic changer by 1 step.
Mag Adv	Advances the magazine by 1 slot.
Mag Rev	Reversely advances the magazine by 1 slot.
Unload Tray	Returns the magazine from the reader to the automatic changer (PULL).
Load Tray	Takes the magazine out of the automatic changer and send it to the reader (PUSH).
START	Starts the measurement.
STOP	Stops the measurement.
Reset Motors	Executes the mechanical initialization operation.

## ■ Menu Bar



It is possible to access various functions of the application software and operate them through the menu bar.

By using the functions of the application, various processes can be executed.

### File

<u>N</u> ew	▶ Batch	Create a new batch	Page 22
<u>E</u> xit UD-7900		Exit from the application	Page 22
Note: This is the only way to end the application.			

### View

Badge Results		“Badge Results” displays all measurement results in the batch.	Page 23
Status Log		Display of unexpected events and error occurrences during processing of the current badge	Page 24
Reports		Tools for data check and analysis will be provided	Page 24
Re-Output Measured Data		Re-output of most recent measurement data output to the external PC	Page 24

### Mode

Reader Check	▶ Loader	Self-diagnosis tool of the loader mechanism	Page 25
	▶ Changer	Self-diagnosis of the automatic changer	Page 25
	▶ Lamp (F3 and F4 Check)	Check tool for Lamp F3 and F4	Page 26

### Parameters

List Parameters		Display/Edit of system parameter	Page 27
Processing Modes		Display/Edit of items relating to measurement process	Page 28
Timing Parameters		Display/Edit of timing parameter	Page 30
Reader Functions		Display/Edit of reader function parameter	Page 32
Calibration Settings		Display/Edit of conversion constant	Page 33
Heat Lamp Power		Display/Control of lamp voltage	Page 34
RS232 Options		Setting of RS232 port	Page 34
Parameter Dump	▶ to Computer	Host computer output of parameter	Page 35
	▶ to File	File output of parameter	Page 35


### Help

<u>A</u> bout UD-7900M		Provision of system information	Page 36
------------------------	--	---------------------------------	---------

## ■ Menu Bar File Command

**New/Batch**

**File/New/Batch Command**



The image shows a dialog box titled "Batch Information". It has three input fields: "Batch ID" (a single-line text box), "Operator" (a single-line text box), and "Comment" (a multi-line text area with a vertical scrollbar on the right). At the bottom of the dialog are two buttons: "Cancel" and "Make Batch".

[File/New/Batch] is the command used to create a new batch and also used to newly measure a TL badge. When [File/New/Batch] is selected from the menu of the main screen, the batch information window will be displayed. When entry for "Batch ID" (8 alphanumeric characters, required), "Operator" and "Comment" is completed and the [Make Batch] button is clicked, a new batch will be created.

After clicking the [Make Batch] button, this batch will automatically be selected as the default and measurement data will be stored in the folder whose name is the name of this batch.

**Exit UD-7900**

**File/Exit UD-7900 Command**

This is the only way to end the application on "File/Exit UD-7900". To prevent from data loss, it will be invalidated during measurement.

## ■ Menu Bar View Commands

### Badge Results

### View/Badge Results Command

[View/Badge Results] displays data of the currently applied batch stored in the “ProcRpt.txt” file. Data items are as shown below and use “Word Pad” to display data.

Item	Description
Title	Process report title
Operational Time:	Date and time when measurement started
Operator:	Batch operator
Comment:	Comment
Header	Data section header
Time	Badge measurement time
Try	Magazine number
SlT	Slot number
BdgNo	Badge number
Typ	Badge type
Err	Badge error code
E1	Element 1 dose value
E2	Element 2 dose value
E3	Element 3 dose value
E4	Element 4 dose value
Drk	Dark count
Ref	Reference element
Corr	Reader sensitivity correction factor
C	Rank C
E	Rank E

### Example of output

ProcRpt.txt														
Batch 20150129 Processing Report														
Operational Time: 01-29-2015 19:07:10 Operator:														
Comment:														
Time	Try	SlT	BdgNo	Typ	Err	E1	E2	E3	E4	Drk	Ref	Corr	C	E
16:00:41	1	1	0001034	7		0.000mSv	0.000mSv	0.000mSv	0.000mSv	4	40	0.955	0	0
16:01:52	1	1	4000053	2		0.782mSv	0.764mSv	0.792mSv	0.854mSv	0	30	0.961	7	0
16:02:29	1	1	4000053	2		0.060mSv	0.046mSv	0.002mSv	0.002mSv	2	30	0.975	7	0
16:02:29	1	2	0000000	0	EE	0.000mSv	0.000mSv	0.000mSv	0.000mSv	0	0	0.975	0	0
16:02:29	1	3	0000000	0	EE	0.000mSv	0.000mSv	0.000mSv	0.000mSv	0	0	0.975	0	0
16:02:51	1	4	4000054	2		1.325mSv	1.094mSv	0.771mSv	0.773mSv	2	24	0.974	7	6
16:03:01	1	5	0000001	2		0.000mSv	0.000mSv	0.000mSv	0.000mSv	0	28	0.949	6	1
16:03:01	1	6	0000000	0	EE	0.000mSv	0.000mSv	0.000mSv	0.000mSv	0	0	0.949	0	0
16:03:23	1	7	4000057	2		1.303mSv	1.275mSv	0.815mSv	0.795mSv	2	30	0.958	7	6
16:03:45	1	8	4000058	2		1.232mSv	1.175mSv	0.758mSv	0.820mSv	0	24	0.953	7	6

**Status Results**[View/Status Results Command](#)

View/ Status Results: Events and errors occurred in the current batch can be checked.

**Example of output**

```
TLDErrorLog.txt
05/05/1998 14:31:45 ** Batch 000000d1 created.
05/05/1998 14:32:21 ** START command issued
05/05/1998 14:33:23 ** STOP command issued
05/05/1998 14:34:00 -> R2 RS232 TX dose data ACK time out
05/05/1998 14:34:03 ** STOP command issued
05/05/1998 14:44:22 ** START command issued
05/05/1998 14:44:22 ** Processing started on batch 000000d1: tray 1.
```

**Reports**[View/Reports Command](#)

Refer to pages 13 - 21 and Appendix-C for details of report.

**Re-Output Measured Data**[View/Re-Output Measured Data Command](#)

Refer to Appendix-D for details of the external PC re-output.




## ■ Menu Bar Mode Commands

### Reader Check/Loader

### Mode/Reader Check/Loader Command


“Reader check/Loader” makes possible to manually control the mechanism section and to check signals. Usually, a repair man uses this for troubleshooting at trouble occurrence.



### Reader Check/Changer

### Mode/Reader Check/Changer Command

“Reader check/Changer” makes possible to manually control the automatic changer and to check signals. Usually, a repair man uses this for troubleshooting at trouble occurrence.



Executes the lamp check and compares the lamp voltage, current and power with reference values.

Badge Read Window

File View Mode Parameters Help

Status Log: Reader Check

- Loader
- Changer
- Lamp(F3 and F4 Check)

03/09/2017 10:51:44 AM START button p  
 03/09/2017 10:51:45 AM Processing start  
 03/09/2017 10:53:08 AM F3 RS232 TX d  
 03/09/2017 10:53:34 AM Processing stop  
 03/09/2017 10:54:26 AM START button pressed.  
 03/09/2017 10:54:26 AM Processing started on batch DefaultBatch: tray 1.  
 03/09/2017 10:59:23 AM STOP button pressed.  
 03/09/2017 10:59:39 AM Processing stopped on batch DefaultBatch.

Batch ID: DefaultBatch Location: 1 - 21 Cor. Factor: 0.237 Lamp Flash Count: 6215 Door:  Open

Location	Badge ID	Type	E/C	Dark Cr	Ref. Ele.	Ref.Light	Element1	Element2	Element3	Element4
1-5	0000305	800	0/0	6	126	0.233	52uSv	26uSv	2uSv	4uSv
1-6	0000306	800	0/0	9	51	0.235	32uSv	57uSv	33uSv	69uSv
1-7	0000307	800	0/0	4	110	0.235	109uSv	51uSv	45uSv	75uSv
1-8	0000308	800	0/0	4	140	0.236	32uSv	76uSv	34uSv	62uSv
1-9	0000309	800	0/0	6	181	0.287	80uSv	39uSv	39uSv	70uSv
1-10	0000310	800	0/0	11	88	0.239	6uSv	19uSv	33uSv	55uSv
1-11	0000311	800	0/0	4	121	0.239	94uSv	56uSv	26uSv	56uSv
1-12	0000312	800	0/0	10	54	0.239	0uSv	19uSv	30uSv	51uSv
1-13	0000313	800	0/0	5	104	0.239	63uSv	69uSv	38uSv	81uSv
1-14	0000314	800	0/0	3	154	0.235	109uSv	89uSv	27uSv	58uSv
1-15	0000315	800	0/0	11	34	0.239	50uSv	31uSv	27uSv	18uSv

Measured Data

Voltage: Lamp V: 14.42 V Lamp B.V: 01.00 V  
 Lamp C: 20.24 A Lamp B.C: 02.28 A  
 Lamp P: 292 W Lamp B.P: 2 W

Mode:  NCR Parity  Anneal Mode

Operation:

Glow Curve

Ready


## ■ Menubar Parameters Commands

### List Parameters


### Parameters/List Parameters Command

“List Parameters” displays and edits parameter P0 - P141. Unused parameters are unavailable to edit. The [Load from File] button can be used to execute backup or load the “parametr.xml” file for work. The [Save to File As] button can be used to execute backup of the current “parametr.xml” file.

The main parameters are grouped and they can be edited on the screens on “Processing Modes”, “Timing Parameters”, “Reader Functions”, “Calibration Settings”, “Heat Lamp Power” and “RS 232 Options”.



Handles all items relating to measurement process together.



**Apply Reader Sensitivity Correction Factor**

Based on reference light read value before correction as reference value, the reader sensitivity correction count will be calculated from the ratio of the value read later on to the reference value. Normally, this is to be applied (having a mark on the check box) for correction of the reader sensitivity. Dirt or dust stuck on the blue filter of PMT may be the possible cause of fluctuation of the reader sensitivity.

**Apply Badge Rank Correction Factors**

TL badge is basically once radiated and measured in the production process at factory and the sensitivity rank is determined based on this measurement result. This rank is given to each TL badge, in form of the rank correction code, by 1 element or by 2 elements and is used to reduce sensitivity fluctuation between element at measurement. Normally, this setting is to be applied (having a mark on the check box).

**Background Dose Subtraction**

For parameter P-4F, set the background to be applied to all elements. The problem when using P-4F is the remaining dose indicated by LiBO element and CaSO element as the result is unequal (LiBO > CaSO). Therefore, an additional background for LiBO can be set for P-45. Generally, the background process is not executed on UD-7900 but on the PC. Thus, it is recommended to set P-45 and P-4F to “0000”.


Table 1 Element definition of TL badge

Dosimeter	Parameter		Hexadecimal		Binary Code for Each Element			
	E1E2	E3E4	E1E2	E3E4	E1	E2	E3	E4
UD-800	13	14	C1C1	E0E0	11000001	11000001	11100000	11100000
UD-801	15	16	C1C1	E0E0	11000001	11000001	11100000	11100000
UD-802	17	18	C1C1	E0E0	11000001	11000001	11100000	11100000
UD-803	19	1A	C1E1	0000	11000001	11100001	00000000	00000000
UD-804	1B	1C	00C0	E0E0	00000000	11000000	11100000	11100000
UD-805	1D	1E	C1E0	0000	11000001	11100000	00000000	00000000
UD-806	1F	20	C1C1	E1E1	11000001	11000001	11100001	11100001
UD-807	21	22	00C1	0000	00000000	11000001	00000000	00000000
UD-808	23	24	C3C3	E0C3	11000011	11000011	11100000	11000011
UD-809	25	26	C3E2	C2E2	11000011	11100010	11100010	11100010
UD-810	27	28	0000	0000	00000000	00000000	00000000	00000000
UD-811	29	2A	0000	C0E0	00000000	00000000	11000000	11100000
UD-812	2B	2C	0000	0000	00000000	00000000	00000000	00000000
UD-813	2D	2E	0000	0000	00000000	00000000	00000000	00000000
UD-814	2F	30	0000	0000	00000000	00000000	00000000	00000000
UD-815	31	32	0000	C0E0	00000000	00000000	11000000	11100000

Table 2 Bit definition of binary code for each element


Bit #	Description
Bit 7	Is this element present in this dosimeter? 0 = no, 1 = yes
Bit 6	Is a rank correction factor to be applied to this element? 0 = no, 1 = yes
Bit 5	Which rank correction factor is to be applied to this element? 0 = C rank, 1 = E rank
Bit 4	Not used (always 0)
Bit 3	Not used (always 0)
Bit 2	Not used (always 0)
Bit 1	Which phosphor sensitivity ratio is used for this element? 00 = CaSO/CaSO (P-3B), 01 = CaSO/nLinBO (P-3C), 10 = CaSO/6Li10BO (P-3D), 11=CaSO/7Li11BO (P-3E)
Bit 0	

All items relating to heating and measurement will be handled together.



Timing parameter (T1-T12)


T1'	Factory setting: LiBO Pre-lamp On time (ms)	Time to replace the heating lamp can be changed
T2'	Factory setting: LiBO Pre-lamp Off time (ms)	Time to replace the heating lamp can be changed
T3'	Factory setting: LiBO Main-lamp On time (ms)	Time to replace the heating lamp can be changed
T0	Factory setting: Pre-heating time (ms)	
T1	Factory setting: CaSO Pre-lamp On time (ms)	Time to replace the heating lamp can be changed
T2	Factory setting: CaSO Pre-lamp Off time (ms)	Time to replace the heating lamp can be changed
T3	Factory setting: CaSO Main-lamp On time (ms)	Time to replace the heating lamp can be changed
T4	Factory setting: Main-lamp On time (ms) - Main count time	
T5	Factory setting: Main-lamp Off time (ms) - Main count time	
T6	Factory setting: Main-lamp Off time (ms)	
T7A	Factory setting: Post-lamp On time (A)	
T8A	Factory setting: Next heating interval time	
T7B	Factory setting: Post-lamp On time (B)	
T8B	Factory setting: Next heating interval time	
T7C	Factory setting: Post-lamp On time (C)	
T8C	Factory setting: Next heating interval time	
T9	Factory setting: Carbon-14 measurement time	
T10	Factory setting: Dark measurement time	
T11	Factory setting: Reference element lamp On time	
T12	Factory setting: Reference element lamp Off time	



## Reader Functions (710 Read Mode)

## Parameters/Reader Functions Command

“Reader Functions (710 Read Mode)” has a collection of frequently-used functions on UD-710 and it is easy to change the reader operation on this screen.



### Check TLD Badge for Parity

Designates presence or absence of parity of punched hole on TL badge. Normally, this setting is to be applied (having a mark on the check box).

### Calculate Reader Sensitivity Value

When this option is designated, the reader will determine the latest reader sensitivity at load of magazine or at start of measurement. At the above timing, Carbon-14 measurement will be executed and the ratio to the reference value will be calculated. This value will be displayed on the screen, stored in P-57 and inserted in transmission data.

### Use Conventional Units of R

When this option is designated, “R” will be used for the unit of dose. When not designated, “Sv” will be used.

### Check Carbon-14 Light Source

When this option is designated, the reader will execute the sensitivity check by Carbon-14 before measuring each TL badge. When the reader sensitivity is out of the reference range, an error occurs and the operation will automatically stop.

### Data String Includes CR and LF

When this option is designated, CR and LF codes will be attached to the end of transmission data.


### Anneal Mode (Use T7C and T8C)

When this option is designated, the anneal time will be extended. “T7C” or “T8C” will automatically be selected according to the measurement results.

## Reader Functions (716 Read Mode)

## Parameters/Reader Functions Command

“Reader Functions (716 Read Mode)” has a collection of frequently-used functions on UD-716 and it is easy to change the reader operation on this screen.



### Basic Mode, Transmit PRE-MAIN-POST “C-string”, Fixed Transmission Data Length

Refer to the “RS232 Options” screen.

### Anneal Mode

When this option is designated, the anneal time will be extended. “T7C” or “T8C” will automatically be selected according to the measurement results.

### Calculate Reader Sensitivity value

When this option is designated, the reader will determine the latest reader sensitivity at load of magazine or at start of measurement. At the above timing, Carbon-14 measurement will be executed and the ratio to the reference value will be calculated. This value will be displayed on the screen, stored in P-57 and inserted in transmission data.

### Use Conventional Units of R



When this option is designated, “R” will be used for the unit of dose. When not designated, “Sv” will be used.

### Subtract Backgrounds, Apply Rank Correction Factor, Apply Reader Sensitivity Factor

Refer to the “Processing Options” screen.

### Check Carbon-14 light Source

When this option is designated, the reader will execute the sensitivity check by Carbon-14 before measuring each TL badge. When the reader sensitivity is out of the reference range, an error occurs and the operation will automatically stop.

Data string include CR and LF

When this option is designated, CR and LF codes will be attached to the end of transmission data.

### External computer control (Expect STX & ETX)

When this option is designated, the format of [STX][ETX][BCC] will be used for reception in the BASIC mode. When not designated, the format will be [CR][LF].

## Calibration Settings

## Parameters/Calibration Settings Command

“Calibration Settings” handles all of each coefficient value together.

Parameter	Value	Unit
[P_3B] CaSO/CaSO	100	
[P_3C] CaSO/nLinBO	1300	
[P_3D] CaSO/6Li10BO	1462	
[P_3E] CaSO/7Li11BO	1400	
[P_3F] CaSO/6LiF	100	
[P_40] CaSO/7LiF	100	
[P_41] CaSO/G	100	
[P_42] CaSO/H	100	
[P_5A] Counter conversion coefficient	500	uSv/(10 <sup>4</sup> count)
[P_57] C-14 Cal. Std.	3200	count
[P_5E] C-14 Avg. 10 meas.	2257	count
[P_71] Reader Sensitivity	0168	
[P_78H] Sensitivity Limit	20	%

### Phosphor Sensitivity Ratios

These parameters are used as a multiplier to the count value of the fluorescent material other than CaSO. The sensitivity ratio of CaSO to each fluorescent material other than CaSO will be displayed.

### Reader Sensitivity

P-57 stores Carbon-14 reference count of the P-counter. Carbon-14 will be read in each time of magazine loading and the read value will be stored in P-67 for the P-counter and in P-68 for the frequency counter. The average value will be stored in P-5E. The sensitivity correction coefficient value of the reader will be calculated by “P-5E/P-57”. When “P-5E/P-57” is within “100%+/- P-78H”, the sensitivity correction coefficient value will be stored in P-71.

The reader will execute calculation of “P-67/P-57” for the sensitivity check. If “P-67/P-57” is within “100% +/- P-78H”, the value will be stored in P-6D.

The sensitivity correction coefficient value will be used for the dose correction. The sensitivity check will be executed in each time of badge measurement, but this value will not be applied to the dose correction.

When cleaned the reader measuring-related section, set a value to P-57 so that P-71 is close to “1”. However, when changed the reference value, it is necessary to execute the correction.

### CCF


The unit used for PCCF is “10-5 mR/P-count”.  
(Do not change it from the factory setting.)

---

<b>Heat Lamp Power</b>	<b>Parameters/Heat Lamp Power Command</b>
------------------------	---

---

“Heat Lamp Power” controls the heating lamp.



### Heat Lamp Power


P-59 stores the reference value of the heating lamp power. If the deviation from the reference value of the actual lamp power (P-6B) exceeded the allowable range (stored in P-77, usually 20%), an error occurrence will be made to interrupt the measurement.

---

<b>RS232 Options</b>	<b>Parameters/RS232 Options Command</b>
----------------------	---

---

The settings relating to the RS232 communication port and the serial printer can be configured.



**RS232 Communication ON**

When this option is selected by putting a mark on the check box, the reader will transmit data to the designated port with the designated conditions (transmission speed, data length, stop bit number, parity, flow control and interface mode).

**RS232 Serial Printer ON**

When this option is selected by putting a mark on the check box, the reader will output data to the printer using the designated port and transmission speed.

**Basic Mode**

When this option is selected by putting a mark on the check box, the reader will be operated in the BASIC mode. In this mode, the reader will request reception check of transmission data from the external device. The external device will transmit ACK if received correctly and transmit NAK to the reader if not.

**Transmit Glow Curve**

When this option is selected by putting a mark on the check box, the reader will output glow data (format code "G").

**Transmit PRE-MAIN-POST**

When this option is selected by putting a mark on the check box, the reader will also output PMP data (format code "C") containing PRE, DOSE-P, DOSE-F and POST in addition to the normal measurement data (format code "G"). This data is useful to determine the crossover point.

**Fixed Transmission Data Length**

When this option is selected by putting a mark on the check box, length of output data in the BASIC mode will be fixed.

**ACK WAIT TIME**

"ACK WAIT TIME" determines the ACK wait time.

**Main Unit COM port**

Set the COM port to be connected with the main unit.

**Parameter Dump to Computer****Parameters/COM Dump to Computer Command**

Sends the parameters to the external device. At that time, the format will be applied by following the interface mode designated in "RS232C Options". Data of this time will be also be stored in the "710pdump.txt" file or the "716pdump.txt" file.

**Parameter Dump to File****Parameters/COM Dump to File Command**


Outputs the parameters as a file.

## ■ Menu Bar Help

**About UD-7900**

**Help/About UD7900 Command**

Displays the application version and the system information.



### 3. How to Operate Windows

---


#### ■ Menu bar navigation

With the menu bar as the starting-point, accessing various functions is available in the common manner of Windows operation.

As an example, see the following procedures of how to display the batch creation screen.


##### Step 1

Search for [File] on the menu bar.



##### Step 2


When the [File] is clicked, the pull-down menu will be displayed.



##### Step 3


Move the mouse pointer onto [New].

“■” indicates that the sub-pull-down menu is available.



##### Step 4

Move the pointer onto “Batch” and click it.



## Appendix-A Error Code

### ■ Error code table

A Automatic Magazine Changer	A3	Error in up/down operation of the magazine shelf
	A5	Operation error to PUSH the magazine to the reader main body or to PULL out of the reader main body
	A6	Solenoid operation error
	A8	Detected existence of magazine in both the measurement section and the magazine changer section
	A9	M1 motor alarm
	AA	M2 motor alarm
	AB	M1 motor (Upper shelf overrun)
	AC	M1 motor (Lower shelf overrun)
B Badge Inserter	B0	Badge push-up operation error
	B1	Badge pull-down operation error
	B2	Magazine return error
	B3	Slider origin error
	B4	Slider end position undetected
	B5	Slider overrun error
	B7	Element position detection error
	BA	1) Element breakage detection error when the card carrier is at the origin point 2) Element breakage detection error
	BE	PM2-related sensor failure
	BF	Irradiation position detection error
C Controller	C1	FPGA start error
	C2	FPGA bus error
	C3	FPGA click cut-off
	C4	System parameter error (checksum error, illegal value)
	C5	System parameter error (Shall be BCD format or 2-byte pack format, but HEX format is applied)
	C7	Power voltage error (15V)
	C8	Front door open
	C9	Power shut-down during element measurement or dose transmission
	CA	PC operation section - Communication error between the control sections
	D Dosimeter	D0
D1		Exceeded the management dose (Element 1, 2)
D2		Exceeded the management dose (All elements)
D4		Badge type out-of-designation
D5		Exceeded the post management dose
H Heating Lamp Driver	H0	Lamp power exceeded the reference value
	H1	Timing parameter error
	H2	Exceeded the lamp flash count (1,000,000 times)

	H3 Error occurred in I2C for halogen lamp voltage current monitoring
I Hybrid Integrator	I2 Out of measurement range at TL badge measurement
	I3 Abnormal value for calibration light source (CAL)
	I4 Exceeded the dark count reference value(both of P-/F-counters)
	I5 Exceeded the reference value
M Magazine Driver	M0 Slot sending error (Cannot send to the next slot)
	M3 Magazine end point error (Magazine is not being sent to the end point position)
N Number Code Reader	N2 NCR parity error
	N3 NCR I2C error
T Temperature	T0 Thermal error around the photomultiplier during measurement
PCI Device Error	Appendix-F Refer to "Measurement and Automation Explorer".

# Appendix-B Example of Report

## Export Data

ExportDose - Notepad

File Edit Format View Help

sCnts	E3MCnts	E3PostCnts	E3Err	E4Dose	E4Cnts	E4Pre	E4Main	E4Post	E4PreCnts	E4MCnts	E4PostCnts	E4Err	Drk	Ele	Light	C	E	Unit
0009	0077	1122	0230	1196	0052	0230	1186	0440	1196	2159	0006	0661	2214	0	2			R
0019	0063	0959	0243	1208	0067	0243	1189	0575	1208	2161	0007	0414	2214	1	1			R
0008	0039	1120	0203	1168	0037	0203	1173	0309	1168	2143	0013	0332	2207	1	1			R
0014	0548	1475	0171	1144	0048	0171	1169	0400	1144	2142	0010	0400	2210	0	1			R
0008	0433	1400	0203	1172	0021	0203	1176	0177	1172	2149	0011	0416	2212	1	1			R
0024	0688	1540	0218	1184	0076	0218	1195	0640	1184	2165	0004	0424	2211	0	1			R
0006	0428	1432	0244	1204	0043	0244	1161	0360	1204	2134	0014	0374	2209	0	1			R
0019	0581	1553	0268	1225	0037	0268	1215	0309	1225	2181	0006	0418	2210	1	1			R
0020	0503	1465	0252	1212	0037	0252	1182	0309	1212	2153	0011	0395	2211	1	1			R
0023	0546	1441	0292	1248	0057	0292	1191	0480	1248	2162	0008	0332	2212	0	0			R
0002	0430	1473	0156	1133	0000	0156	1195	0000	1133	2165	0018	0379	2212	1	1			R
0031	0593	1530	0316	1265	0037	0316	1170	0309	1265	2143	0008	0391	2210	1	2			R
0038	0584	1454	0182	1155	0042	0182	1182	0354	1155	2155	0011	0381	2213	1	2			R
0021	0507	1373	0296	1248	0053	0296	1226	0442	1248	2189	0009	0352	2209	1	2			R
0010	0425	1425	0256	1212	0043	0256	1201	0360	1212	2167	0007	0392	2207	0	1			R
0024	0534	1548	0319	1270	0099	0319	1214	0840	1270	2181	0005	0379	2212	1	1			R
0007	0528	1500	0203	1168	0048	0203	1195	0400	1168	2161	0009	0364	2207	0	0			R
0023	0454	1427	0378	1314	0069	0378	1217	0575	1314	2180	0008	0398	2208	1	1			R
0022	0394	1465	0167	1141	0037	0167	1149	0309	1141	2126	0009	0298	2212	1	1			R
0020	0475	1458	0190	1160	0033	0190	1167	0280	1160	2140	0007	0352	2210	0	2			R
0034	0469	1453	0187	1156	0053	0187	1151	0440	1156	2126	0006	0335	2209	0	1			R
0017	0448	1454	0280	1234	0058	0280	1215	0486	1234	2179	0006	0335	2209	1	1			R
0011	0441	1492	0286	1236	0073	0286	1200	0600	1236	2165	0007	0401	2207	0	1			R
0007	0506	1417	0282	1236	0053	0282	1152	0440	1236	2128	0012	0360	2210	0	2			R
0023	0473	1443	0215	1181	0052	0215	1166	0434	1181	2139	0007	0384	2210	7	1			R

## Badge Results

ProcRpt - Notepad

File Edit Format View Help

Batch 500Badge Processing Report

Operational Time: 12/13/2016 4:40:19 PM Operator: test  
Comment: test

Time	Try	Slit	BadgeNo	Typ	Err	E1	E2	E3	E4	Drk	Ref	Corr	C	E
16:40:20	01	01	00000000	08		122uSv	169uSv	9uSv	230uSv	006	661	0.207	0	2
16:40:53	01	02	00000000	08		145uSv	135uSv	7uSv	243uSv	007	414	0.207	1	1
16:41:13	01	03	00000000	08		101uSv	155uSv	5uSv	203uSv	013	332	0.207	1	1
16:41:33	01	04	00000000	08		157uSv	162uSv	65uSv	171uSv	010	400	0.207	0	1
16:41:54	01	05	00000000	08		89uSv	151uSv	51uSv	203uSv	011	416	0.207	1	1
16:42:15	01	06	00000000	08		171uSv	175uSv	82uSv	218uSv	004	424	0.207	0	1
16:42:35	01	07	00000000	08		91uSv	91uSv	51uSv	244uSv	014	374	0.207	0	1
16:42:56	01	08	00000000	08		105uSv	179uSv	69uSv	268uSv	006	418	0.207	1	1
16:43:16	01	09	00000000	08		136uSv	131uSv	60uSv	252uSv	011	395	0.207	1	1
16:43:36	01	10	00000000	08		179uSv	118uSv	64uSv	292uSv	008	332	0.207	0	0
16:43:56	01	11	00000000	08		89uSv	73uSv	51uSv	156uSv	018	379	0.207	1	1
16:44:16	01	12	00000000	08		185uSv	111uSv	71uSv	316uSv	008	391	0.207	1	2
16:44:36	01	13	00000000	08		99uSv	130uSv	69uSv	182uSv	011	381	0.207	1	2
16:44:56	01	14	00000000	08		174uSv	85uSv	61uSv	296uSv	009	352	0.207	1	2
16:45:15	01	15	00000000	08		155uSv	193uSv	51uSv	256uSv	007	392	0.207	0	1
16:45:36	01	16	00000000	08		146uSv	115uSv	63uSv	319uSv	005	379	0.207	1	1
16:45:57	01	17	00000000	08		121uSv	92uSv	64uSv	203uSv	009	364	0.207	0	0
16:46:18	01	18	00000000	08		101uSv	160uSv	55uSv	378uSv	008	398	0.207	1	1
16:46:37	01	19	00000000	08		152uSv	120uSv	47uSv	167uSv	009	298	0.207	1	1
16:46:57	01	20	00000000	08		110uSv	138uSv	57uSv	190uSv	007	352	0.207	0	2



# Tray Summary

Tray01Rpt - Notepad

File Edit Format View Help

Batch: DefaultBatch2 Tray Report  
 Operational Time: 12/15/2016 5:31:59 PM Operator:  
 Comment:

Sl#	BadgeNo	E1	E2	E3	E4	Average	Sigma
01	0000000	122uSv	169uSv	9uSv	230uSv	132uSv	80.93
02	0000000	145uSv	135uSv	7uSv	243uSv	132uSv	83.85
03	0000000	101uSv	155uSv	5uSv	203uSv	116uSv	73.55
04	0000000	157uSv	162uSv	65uSv	171uSv	138uSv	42.87
05	0000000	89uSv	151uSv	51uSv	203uSv	123uSv	58.14
06	0000000	171uSv	175uSv	82uSv	218uSv	161uSv	49.46
07	0000000	91uSv	91uSv	51uSv	244uSv	119uSv	73.85
08	0000000	105uSv	179uSv	69uSv	268uSv	155uSv	76.22
09	0000000	136uSv	131uSv	60uSv	252uSv	144uSv	68.83
10	0000000	179uSv	118uSv	64uSv	292uSv	163uSv	84.74
...							
40	0000000	165uSv	80uSv	52uSv	293uSv	147uSv	93.75
41	0000000	137uSv	158uSv	59uSv	274uSv	157uSv	76.96
42	0000000	191uSv	187uSv	55uSv	292uSv	181uSv	84.16
43	0000000	158uSv	138uSv	68uSv	224uSv	147uSv	55.61
44	0000000	210uSv	167uSv	65uSv	247uSv	172uSv	68.09
45	0000000	131uSv	113uSv	57uSv	179uSv	120uSv	43.65
46	0000000	127uSv	137uSv	65uSv	245uSv	143uSv	64.77
47	0000000	155uSv	165uSv	41uSv	257uSv	154uSv	76.65
48	0000000	207uSv	185uSv	65uSv	344uSv	200uSv	99.04
49	0000000	126uSv	101uSv	55uSv	266uSv	137uSv	78.71
50	0000000	150uSv	107uSv	59uSv	238uSv	138uSv	65.85
-----							
No. of Data	050	050	050	050			
Average	137uSv	131uSv	54uSv	246uSv			
Sigma	31.21	33.66	14.29	47.64			

# Status

StatusLog - Notepad

File Edit Format View Help

```

2016/12/13 16:40:19 START button pressed.
2016/12/13 16:40:20 Processing started on batch 500Badge: tray 1.
2016/12/13 16:57:53 Processing started on batch 500Badge: tray 2.
2016/12/13 17:15:36 Processing started on batch 500Badge: tray 3.
2016/12/13 17:33:14 Processing started on batch 500Badge: tray 4.
2016/12/13 17:50:57 Processing started on batch 500Badge: tray 5.
2016/12/13 18:08:34 Processing started on batch 500Badge: tray 6.
2016/12/13 18:26:37 Processing started on batch 500Badge: tray 7.
2016/12/13 18:44:41 Processing started on batch 500Badge: tray 8.
2016/12/13 19:02:49 Processing started on batch 500Badge: tray 9.
2016/12/13 19:20:35 Processing started on batch 500Badge: tray 10.
2016/12/13 19:38:46 Processing stopped on batch 500Badge.
2016/12/13 19:48:48 UD-7900 Program Shutdown.
  
```

# Appendix-C Communication Specifications

---

## Basic specifications

The following are descriptions about the communication specifications between UD-7900 and an external device (host computer, printer, etc.).

RS-232C is used for communication.

Communication: RS-232C

Communication speed: 1200-115200 baud rate

Communication method: Half-duplex

Synchronization: Start-stop synchronization

Character code: ASCII 7 bits + parity

Time-out: 3 sec (Changeable by parameter)


Error detection: Parity (Character unit), BCC (Block unit)

Error correction: Re-send

Connection type: 1:1

Communication mode: Basic mode/No Data Link mode

### Character configuration:



### Bit weight:

Bit	b1 b2 b3 b4 b5 b6 b7
Weight	(LSB) 1 2 4 8 16 32 64 (MSB)

### Telegram configuration:

- (1) Control sequence – ACK(06), NAK(21), DC1(17), DC2(18), DC3(19), DC4(20)
- (2) Data telegram

Pause	STX(02)	Data section	ETX(03)	BCC	Pause
-------	---------	--------------	---------	-----	-------

## 710/791 Interface mode

TLD measurement data (Format code '2')

0	stx
1	2
2	number of characters
3	
4	reader code
5	(device address)
6	text number (00-99)
7	
8	space
9	magazine slot number
10	
11	badge code type
12	badge code e-rank
13	badge codec-rank
14	
15	
16	
17	personnel id code
18	
19	
20	
21	element 1 measured data
22	
23	
24	
25	element 2 measured data
26	
27	
28	
29	element 3 measured data
30	
31	
32	
33	element4 measured data
34	
35	
36	
37	dark counts
38	
39	
40	reference element counts
41	
42	
43	reader sensitivity
44	
45	
46	
47	space
48	error code (1)
49	
50	error code (2)
51	
52	error code (3)
53	
54	error code (4)
55	
56	ETX
57	BCC

CAL data (Format code '1')

0	stx
1	1
2	number of characters
3	
4	reader code
5	(device address)
6	text number (00-99)
7	
8	S
9	status
10	
11	magazine number (01-10)
12	
13	Wear-term year
14	
15	wear-term month
16	reference level for element 1
17	(control dose)
18	
19	
20	reference level for element 2
21	(control dose)
22	
23	
24	background level by natural dose
25	
26	
27	
28	sensitivity correction factor
29	
30	
31	
32	operating hour (0000-9999)
33	
34	
35	
36	operating minute (00-59)
37	
38	lamp flash times
39	
40	
41	
42	
43	
44	
45	
46	operating year (lower 2 digit)
47	
48	Operating month (01-12)
49	
50	operating date (01-31)
51	
52	space
53	error code (1)
54	
55	error code (2)
56	
57	error code (3)
58	
59	error code (4)
60	
61	etx
62	bcc



## 716 Interface mode

TLD measurement data (Format code '2')

0	stx
1	2
2	number of characters
3	
4	reader code
5	(device address)
6	text number (00-99)
7	
8	Space
9	magazine slot number
10	
11	badge code type
12	badge code e-rank
13	badge codec-rank
14	personnel id code
15	
16	
17	
18	
19	
20	
21	element 1 measured data
22	
23	
24	
25	element 2 measured data
26	
27	
28	
29	element 3 measured data
30	
31	
32	
33	element4 measured data
34	
35	
36	
37	dark counts
38	
39	
40	reference element counts
41	
42	
43	reader sensitivity check factor
44	
45	
46	
47	p-f flag
48	'Y' (716 mark)
49	sequence number
50	
51	
52	error code (1)
53	
54	error code (2)
55	
56	error code (3)
57	
58	error code (4)
59	
60	etx
61	bcc

CAL data (Format code '1')

0	stx
1	1
2	number of characters
3	
4	reader code
5	(device address)
6	text number (00-99)
7	
8	Space
9	measuring year
10	
11	measuring month
12	
13	measuring day
14	
15	space
16	
17	reference light counts
18	
19	
20	
21	background level by natural dose
22	
23	
24	
25	sensitivity correction
26	factor
27	
28	
29	high voltage
30	
31	
32	lamp voltage
33	
34	
35	lamp flash times
36	
37	
38	
39	
40	
41	
42	
43	space
44	
45	reader function
46	
47	
48	
49	'Y' (716 mark)
50	p-counter conversion coeff
51	
52	
53	
54	f-counter conversion coeff
55	
56	
57	
58	error code (1)
59	
60	error code (2)
61	
62	error code (3)
63	
64	error code (4)
65	
66	etx
67	bcc



## Appendix-D Parameter list

Adr.	Item	Format	Default	Description
00	T0	BCD	0001	Lamp on time(ms).
01	T1	BCD	0065	Lamp on time(ms) for Pre.
02	T2	BCD	0500	Lamp off time(ms) for Pre.
03	T3	BCD	0050	Lamp on time(ms) for Pre.
04	T4	BCD	0045	Lamp on time(ms) for Main.
05	T5	BCD	0400	Lamp off time(ms) for Main.
06	T6	BCD	0150	Lamp off time(ms) for Post.
07	T7A	BCD	0035	Lamp on time(ms) for Post. DOSE < CMPDS1
08	T8A	BCD	1700	Lamp off time(ms) for Post. DOSE < CMPDS1
09	T7B	BCD	0045	Lamp on time(ms) for Post. CMPDS1 ≤ DOSE < CMPDS2
0A	T8B	BCD	1700	Lamp off time(ms) for Post. CMPDS1 ≤ DOSE < CMPDS2
0B	T7C	BCD	0055	Lamp on time(ms) for Post. CMPDS2 ≤ DOSE
0C	T8C	BCD	1700	Lamp off time(ms) for Post. CMPDS2 ≤ DOSE
0D	T9	BCD	1000	Reference Right count time.
0E	T10	BCD	0445	Dark count time.
0F	T11	BCD	0240	Lamp on time(ms) for Reference Element.
10	T12	BCD	0600	Lamp off time(ms) for Reference Element.
11	DARKLC	BCD	0225	Coefficient for converting to DARK value for CAL 1 time (T9). Fixed point (99.99).
12	DARKEC	BCD	0189	A factor for calculating the DARK value for the Reference Element. Fixed point (99.99).
13.H	TYP800E1	HEX	C1	<p>Element definition.</p> <p>The meaning of each bit of element information is as follows. (i = element numbers 1 to 4)</p> <p>Ei [b7] : Presence / absence of element 0 = No, 1 = Yes</p> <p>Ri [b6] : With or without rank correction 0 = No, 1 = Yes</p> <p>Ci [b5] : Selection of rank code 0=upper row (C) rank, 1=lower row (E) rank</p> <p>li [b4] : Presence or absence of I rank correction. This bit is effective only when P51.07 = 1. Function that has rank code (7th digit &amp; 0x7 of badge ID) instead of Ci 0 = No, 1 = Yes</p> <p>[b3] : Reserve</p> <p>CEFn [b2 to 0] : Conversion constant number (CEF0 to 7) 0 to 7 = CEF0 to ECF7</p>
13.L	TYP800E2	HEX	C1	
14.H	TYP800E3	HEX	E0	
14.L	TYP800E4	HEX	E0	
15.H	TYP801E1	HEX	C1	
15.L	TYP801E2	HEX	C1	
16.H	TYP801E3	HEX	E0	
16.L	TYP801E4	HEX	E0	
17.H	TYP802E1	HEX	C1	
17.L	TYP802E2	HEX	C1	
18.H	TYP802E3	HEX	E0	
18.L	TYP802E4	HEX	E0	
19.H	TYP803E1	HEX	C1	
19.L	TYP803E2	HEX	E1	
1A.H	TYP803E3	HEX	00	
1A.L	TYP803E4	HEX	00	
1B.H	TYP804E1	HEX	00	
1B.L	TYP804E2	HEX	C0	
1C.H	TYP804E3	HEX	E0	
1C.L	TYP804E4	HEX	E0	
1D.H	TYP805E1	HEX	00	
1D.L	TYP805E2	HEX	C0	
1E.H	TYP805E3	HEX	E0	
1E.L	TYP805E4	HEX	E0	
1F.H	TYP806E1	HEX	C1	
1F.L	TYP806E2	HEX	E0	

20.H	TYP806E3	HEX	00		
20.L	TYP806E4	HEX	00		
21.H	TYP807E1	HEX	C1		
21.L	TYP807E2	HEX	C1		
22.H	TYP807E3	HEX	E1		
22.L	TYP807E4	HEX	E1		
23.H	TYP808E1	HEX	00		
23.L	TYP808E2	HEX	C1		
24.H	TYP808E3	HEX	00		
24.L	TYP808E4	HEX	00		
25.H	TYP809E1	HEX	C3		
25.L	TYP809E2	HEX	C3		
26.H	TYP809E3	HEX	E0		
26.L	TYP809E4	HEX	C3		
27.H	TYP810E1	HEX	C3		
27.L	TYP810E2	HEX	E2		
28.H	TYP810E3	HEX	E2		
28.L	TYP810E4	HEX	E2		
29.H	TYP811E1	HEX	C3		
29.L	TYP811E2	HEX	C3		
2A.H	TYP811E3	HEX	E2		
2A.L	TYP811E4	HEX	80		
2B.H	TYP812E1	HEX	00		
2B.L	TYP812E2	HEX	00		
2C.H	TYP812E3	HEX	C0		
2C.L	TYP812E4	HEX	E0		
2D.H	TYP813E1	HEX	00		
2D.L	TYP813E2	HEX	00		
2E.H	TYP813E3	HEX	00		
2E.L	TYP813E4	HEX	00		
2F.H	TYP814E1	HEX	00		
2F.L	TYP814E2	HEX	00		
30.H	TYP814E3	HEX	00		
30.L	TYP814E4	HEX	00		
31.H	TYP815E1	HEX	00		
31.L	TYP815E2	HEX	00		
32.H	TYP815E3	HEX	C0		
32.L	TYP815E4	HEX	E0		
33	RNK0	BCD	1000	Rank correction factor 0.	
34	RNK1	BCD	1105	Rank correction factor 1.	
35	RNK2	BCD	1222	Rank correction factor 2.	
36	RNK3	BCD	1350	Rank correction factor 3.	
37	RNK4	BCD	1000	Rank correction factor 4.	
38	RNK5	BCD	0741	Rank correction factor 5.	
39	RNK6	BCD	0819	Rank correction factor 6.	
3A	RNK7	BCD	0905	Rank correction factor 7.	
3B	CEF0	2byte pack	0100	CaSO4 :Tm 0100 (fixed)	Sensitivity coefficient A value 100 times the reciprocal of the sensitivity to CaSO4 is set as the
3C	CEF1	2byte pack	1450	Li2B4O7 :Cu	
3D	CEF2	2byte pack	1692	6Li2B4O7 CEF1/0.65	
3E	CEF3	2byte pack	1600	Li2B4O7:Cu CEF1/0.75	



3F	CEF4	2byte pack	0100	E Reserve	coefficient of each material. That is, the relative ratio when CaSO <sub>4</sub> is set to 100.
40	CEF5	2byte pack	0100	F Reserve	
41	CEF6	2byte pack	0100	G Reserve	
42	CEF7	2byte pack	0100	H Reserve	
43	CSOGND	2byte pack	0000	Background dose of CaSO. (unit of measurement is applied for dose unit)	
44	-	-	-	-	
45	LBOGND	2byte pack	0000	Background dose of LiBO. (unit of measurement is applied for dose unit)	
46	T1_LI	BCD	0065	T1 for LiBO.	
47	T2_LI	BCD	0500	T2 for LiBO.	
48	T3_LI	BCD	0050	T3 for LiBO.	
49.H	YOBI	HEX	-	Reserve	
49.L	MRYCKB	HEX	-	Memory check code.	
4A					
4B	CMPDS1	2byte pack	1100	Post annealing switching value 1, 2. (Reference value for selecting T7, T8) (unit of measurement is applied for dose unit)	
4C	CMPDS2	2byte pack	1500		
4D	-	-	-	-	
4E	-	-	-	-	
4F	BKGND	2byte pack	0000	Background dose. (unit of measurement is applied for dose unit)	
50.15	TXDATEDATA	bool	false	Add time information to measurement data (Format 2) (true) or not (false).	
50.14	BASIC	bool	false	BASIC mode (true), No data link mode (false).	
50.13	-	-	-	-	
50.12	ANNEALMODE	bool	false	Annealing measurement mode (true), normal measurement mode (false).	
50.11	NCRPT	bool	true	Check the parity of NCR (true), do not (false)	
50.10	CALCC14	bool	true	Updates the reader sensitivity correction coefficient[P71] during CAL measurement (true) and does not (false).	
50.09	TXPMPDATA	bool	false	Outputs PRE/MAIN/POST to RS232 (true) or not (false).	
50.08	-	-	-	-	
50.07	USEREM	bool	false	The dose is R unit (true), Sv unit (false).	
50.06	-	-	-	-	
50.05	-	-	-	-	
50.04	-	-	-	-	
50.03	BGCOR	bool	false	Background correction (true), not (false).	
50.02	BDCOR	bool	true	Apply the rank correction coefficient (true), do not (false).	
50.01	-	-	-	-	
50.00	RDCOR	bool	true	Apply the reader sensitivity correction factor [P71] (true), do not (false).	
51.15	-	-	-	-	
51.14	STX_ETX	bool	false	Set the reception in BASIC mode to [STX] [ETX] [BCC] format (true), set to [CR] [LF] format (false).	
51.13	-	-	-	-	
51.12	-	-	-	-	
51.11	-	-	-	-	
51.10	-	-	-	-	
51.09	-	-	-	-	
51.08	-	-	-	-	
51.07	RNK_I	bool	false	Rank correction follows (true), not obey (false) according to li	

				[b4] of Badge Table [P13 to 32].
51.06	CHECKC14	bool	true	CAL operation switching. (True): [START] and 10 times when magazine is inserted (store the average in [P5E]), once every badge measurement (SNS unsaved). (False): [START] and 5 times when magazine is inserted (store accumulation in [P5E]), 0 times for badge measurement (without CAL).
51.05	-	-	-	-
51.04	-	-	-	-
51.03	-	-	-	-
51.02	ADD_CR	bool	false	[CR] [LF] is added (true) or not (false) to the transmission message in BASIC mode.
51.01	-	-	-	-
51.00	DATA_FIXED	bool	false	Send transmission message in BASIC mode with fixed length (64 bytes) (true) or not (false).
52	-	-	-	-
53	SREFEL	2byte pack	0300	Ref. Element reference value (count). $(P53 - P56) \leq$ Ref. Ele. $<(P53 + P56)$ I5 error if it is out of range.
54	SDARKP	2byte pack	0300	DARK reference value (count). I4 error when exceeding.
55	-	-	-	-
56	SREFER	2byte pack	0300	Ref. Element Check width (count). $(P53 - P56) \leq$ Ref. Ele. $<(P53 + P56)$ I5 error if it is out of range.
57	SNSSTD	2byte pack	2500	Ref. Light (CAL) reference value (count). $P57 - P57 \times P78.H \div 100 \leq P67$ (CAL value) $\leq P57 + P57 \times P78.H \div 100$
58	-	-	-	-
59	LPWST	BCD	0300	Heating lamp power reference value. If it is outside the range, H0 error.
5A	CNVCEP	BCD	0160	Conversion constant. Dose per 4th power of 10 (uSv). In the case of R units, it is interpreted as the dose per 10 5 power (mR).
5B	-	-	-	-
5C	-	-	-	-
5D	-	-	-	-
5E	SNS	2byte pack	-	Ref. Light (CAL) reading (count). When P51.06 = 1 Update with [CALA] 10 times CAL average value at [START] and magazine insertion. However, it does not update when a CAL error occurs. When P51.06 = 0 It is updated with [CAL integrated value 5 times at the time of [START] and magazine insertion. However, it does not update when a CAL error occurs.
5F	DARKP	2byte pack	-	DARK reading (count). DARK count value for T10. I4 error when exceeding.
60	-	-	-	-
61	DARKEP	2byte pack	-	DARK value (for Ref. Ele.) (Count). = $P5F \times (T11 + T12) \div T10$
62	-	-	-	-
63	REFELP	2byte pack	-	Ref.Element value (count). (T11 + T12)
64	-	-	-	-
65	DARKLP	2byte pack	-	DARK value (for Ref. Light) (count). = $P5F \times T9 \div T10$
66	-	-	-	-

67	REFLP	2byte pack	—	Ref. Light value (count). (T9)
68	-	-	-	-
69	-	-	-	-
6A	LMPVLT	BCD	—	Heating lamp voltage reading. Fixed point (99.99 V)
6B	LMPPWR	BCD	—	Heating lamp power value. H0 error if it is out of range.
6C	LMPCUR	BCD	—	Heating lamp current reading. Fixed point (99.99 A)
6D	B_LMPVLT	BCD	—	Heating lamp voltage reading at lamp OFF. Fixed point (99.99 V)
6E	B_LMPPWR	BCD	—	Calculated heating lamp power at lamp OFF.
6F	B_LMPCUR	BCD	—	Heating lamp current reading value when lamp is off. Fixed point (99.99 A)
70	-	-	-	-
71	SCEF	BCD	1000	Leader sensitivity correction coefficient. Fixed decimal point (9.999)
72.H	ELEMENT	HEX	0F	Designated element (00 to 0F). ELEMENT: Necessity of measurement for each element number. When 0 is set to the relevant bit, measurement of the corresponding element number is prohibited. Applicable Bit (arrangement) → Bit 3 to 0 = Elements 1 to 4
72.L	-	-	-	-
73	-	-	-	-
74	-	-	-	-
75	-	-	-	-
76.H	ACKTIM	BCD	03	ACK wait time. (sec)
76.L	TRTMR	BCD	10	Transmission delay timer. (100 msec unit)
77.H	-	-	-	-
77.L	LPWER	BCD	20	Heating lamp power check width(%). If it is outside the range, H0 error.
78.H	SNSER	BCD	20	Ref. Light check width(%). I3 error if it is out of range.
78.L	RDERCD	BCD	01	Reader number.
79	-	-	-	-
7A	-	-	-	-
7B	-	-	-	-
7C	FLASH	BCD	00000000	Number of times the heating lamp was flashed. H2 error when exceeding.
7D				
7E.H	PROM_NO	BCD	—	Program No. Unused.
7E.L	PRONO	BCD	—	Production number (6 digits). Unused.
7F.H				
80	-	-	-	-
81	-	-	-	-
82	-	-	-	-
83	FPGAVER	BCD	—	FPGA program Ver. Unused.
84	CPLDVER	BCD	—	CPLD program Ver. Unused.
85	-	-	-	-
86	-	-	-	-
87	-	-	-	-
88	-	-	-	-
89	-	-	-	-
8A	-	-	-	-
8B	-	-	-	-
8C	-	-	-	-

8D	-	-	-	-
8E	-	-	-	-
8F	-	-	-	-
90	RFSKIN	2byte pack	1500	E1 controlled dose. D1 error when exceeding.(unit of dosage applied at the time of measurement)
91	RFBODY	2byte pack	1500	E2 controlled dose. D1 error when exceeding.(unit of dosage applied at the time of measurement)
92	RFANL	2byte pack	1500	Controlled dose of all elements. D2 error when exceeding. (Unit of measurement is applied for dose unit)
93	RFPOST	2byte pack	1100	Controlled Post dose of all elements. D5 error when exceeding. (Unit of measurement is applied for dose unit)
94	RLCEF	BCD	—	Badge sensitivity. (= Ref. Light 1 - Readout ÷ Standard) Fixed point (9.999).
95	-	-	-	-
96.H	TRAY	HEX	—	TRAY (Magazine) No.
96.L	SLOT	HEX	—	Slot No.
97	-	-	-	-
98	-	-	-	-
99	BADGE	HEX	7FFF	Badge Assignment (0000 - FFFF). D4 error if unspecified Badge. Specify validity / invalidity of UD-800 (b 0) to 815 (b 15) with 16 bits.
9A	OVRDS1	2byte pack	0000	LiBO limit value.
9B	OVRDS2	2byte pack	0000	CaSO limit value.
9C	-	-	-	-
9D	-	-	-	-
9E	-	-	-	-
9F	-	-	-	-
A0.15	DCCOR	bool	true	Dark Correct (true) or not (false).
A0.14	-	-	-	-
A0.13	-	-	-	-
A0.12	ELEMENTCH	bool	true	Check element break (true), not (false). D0 error when element breaks.
A0.11	CYCLELOC	bool	true	Restore Location at startup (true), not (false).
A0.10	DOORCHK	bool	true	Door check (true), do not (false). Unused.
A0.09	-	-	-	-
A0.08	-	-	-	-
A0.07	-	-	-	-
A0.06	-	-	-	-
A0.05	-	-	-	-
A0.04	-	-	-	-
A0.03	-	-	-	-
A0.02	-	-	-	-
A0.01	-	-	-	-
A0.00	-	-	-	-
A1	-	-	-	-
A2	-	-	-	-
A3	-	-	-	-
A4	-	-	-	-
A5	-	-	-	-
A6	-	-	-	-
A7	-	-	-	-
A8	-	-	-	-

A9	-	-	-	-
AA	-	-	-	-
AB	-	-	-	-
AC	-	-	-	-
AD.H	ACKTRN	BCD	02	ACK retry count. If retry over, R1 error.
AD.L	-	-	-	-
AE	-	-	-	-
AF	-	-	-	-
B0	-	-	-	-
B1	-	-	-	-
B2	-	-	-	-
B3	-	-	-	-
B4	-	-	-	-
B5	-	-	-	-
B6	-	-	-	-
B7	PMP1PRE	2byte pack	—	E1 Pre dose.
B8	PMP1MAIN	2byte pack	—	E1 Main dose.
B9	PMP1POST	2byte pack	—	E1 Post dose.
BA	PMP2PRE	2byte pack	—	E2 Pre dose.
BB	PMP2MAIN	2byte pack	—	E2 Main dose.
BC	PMP2POST	2byte pack	—	E2 Post dose.
BD	PMP3PRE	2byte pack	—	E3 Pre dose.
BE	PMP3MAIN	2byte pack	—	E3 Main dose.
BF	PMP3POST	2byte pack	—	E3 Post dose.
C0	PMP4PRE	2byte pack	—	E4 Pre dose.
C1	PMP4MAIN	2byte pack	—	E4 Main dose.
C2	PMP4POST	2byte pack	—	E4 Post dose.
C3	-	-	-	-
C4	-	-	-	-
C5	FULPMP1PREP	2byte pack	—	E1 Pre count.
C6	-	-	-	-
C7	FULPMP1MAINP	2byte pack	—	E1 Main count.
C8	-	-	-	-
C9	FULPMP1POSTP	2byte pack	—	E1 Post count.
CA	-	-	-	-
CB	FULPMP2PREP	2byte pack	—	E2 Pre count.
CC	-	-	-	-
CD	FULPMP2MAINP	2byte pack	—	E2 Main count.
CE	-	-	-	-
CF	FULPMP2POSTP	2byte pack	—	E2 Post count.
D0	-	-	-	-
D1	FULPMP3PREP	2byte pack	—	E3 Pre count.
D2	-	-	-	-
D3	FULPMP3MAINP	2byte pack	—	E3 Main count.
D4	-	-	-	-
D5	FULPMP3POSTP	2byte pack	—	E3 Post count.
D6	-	-	-	-
D7	FULPMP4PREP	2byte pack	—	E4 Pre count.
D8	-	-	-	-

D9	FULPMP4MAINP	2byte pack	—	E4 Main count.
DA	-	-	-	-
DB	FULPMP4POSTP	2byte pack	—	E4 Post count.
DC	-	-	-	-
DD	-	-	-	-
DE	-	-	-	-
DF	-	-	-	-
E0	DEADTIME	BCD	0200	Dead time of photomultiplier tube (0.1 ns unit).
E1	-	-	-	-
E2	-	-	-	-
E3	-	-	-	-
E4	PROCOT	HEX	—	Count value.
E5				
E6	-	-	-	-
E7	-	-	-	-
E8	-	-	-	-
E9	-	-	-	-
EA	-	-	-	-
EB	-	-	-	-
EC	-	-	-	-
ED	-	-	-	-
EE	-	-	-	-
EF	-	-	-	-
F0	-	-	-	-
F1	-	-	-	-
F2	-	-	-	-
F3.15	-	-	-	-
F3.14	-	-	-	-
F3.13	-	-	-	-
F3.12	-	-	-	-
F3.11	-	-	-	-
F3.10	-	-	-	-
F3.09	-	-	-	-
F3.08	-	-	-	-
F3.L	-	-	-	-
F4	-	-	-	-
F5	-	-	-	-
F6	SENDGLOWTHR ESHOLD	2byte pack	0000	Glow data output threshold. (Dose unit applies unit of measurement point)
F7	-	-	-	-
F8	-	-	-	-
F9	-	-	-	-
FA.15	-	-	-	-
FA.14	-	-	-	-
FA.13	-	-	-	-
FA.12	-	-	-	-
FA.11	-	-	-	-
FA.10	-	-	-	-
FA.09	-	-	-	-
FA.08	TXGLOWDATA	bool	false	Output glow data (true), do not output (false).
FA.L	-	-	-	-

FB	-	-	-	-
FC	-	-	-	-
FD	-	-	-	-
FE	-	-	-	-
FF	-	-	-	-

**Panasonic System Solutions Japan Co., Ltd.**

600, Saedo-cho, Tsuzuki-ku, Yokohama-shi, Kanagawa 224-8539, Japan

PYQX1101ZA/J1  
Jan. 2017  
Published in Japan